

BODY ASSIST 988 THERMAL ELBOW WRAP



- One size thermic elbow wrap for simple elbow related strains and sprains.
- Medical grade CR neoprene is LATEX FREE.
- Medical plush inner liner draws perspiration, delivers warmth and provides soft plush comfort against the skin.
- Superior compression.
- Adjustable closure with contoured design for a snug fit that easily fits under most shirts.
- Black looks great as part of any sports uniform.
- Suited to minor sports injuries.
- Fits left or right.
- Washable & hygienic.
- Simple one hand application.
- One size fit. Application:
- To provide assisted relief of tendonitis, bursitis, sprains, strains and arthritis.
- Retains body warmth & can increase blood circulation.

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name
FBS307770	BODY ASSIST 988 THERMAL ELBOW WRAP - UNIVERSAL