

## Venosan - 'Microfibre' Compression Stockings



- Microfibre is designed to provide relief from Mild Oedema, Slight Varicose Veins and Tired Legs
- **Styles:** Knee High. All stockings are closed-toe only
- **Colours:** Male: Black and Female: Black and Beige
- **Compression:** Preventative to Low-Therapeutic range: 15 - 20mmHg (male & female) and 18 - 22mmHG (male only)

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Supporting and stimulating blood circulation

**Category:** Low risk

Click [here](#) for more information

Code	Display Name	Colour
PTS641965	Venosan 5001 (Cotton) - Knee - 18-21mmHg - Unisex - O/Toe - Sml - Marokko	Marokko
PTS641996	Venosan - Microfibre - Knee - 15-20mmHg - Male - C/Toe - Sml - Blk	Black
PTS641965BLL	Venosan 5001 (Cotton) - Knee - 18-21mmHg - Unisex - C/Toe - Lge - Blk	Black
PTS641996L	Venosan - Microfibre - Knee - 15-20mmHg - Male - C/Toe - Lge - Blk	Black
PTS641965BLM	Venosan 5001 (Cotton) - Knee - 18-21mmHg - Unisex - C/Toe - Med - Blk	Black
PTS641996M	Venosan - Microfibre - Knee - 15-20mmHg - Male - C/Toe - Med - Blk	Black
PTS641965BLS	Venosan 5001 (Cotton) - Knee - 18-21mmHg - Unisex - C/Toe - Sml - Blk	Black

Code	Display Name	Colour
PTS641996XL	Venosan - Microfibre - Knee - 15-20mmHg - Male - C/Toe - XL - Blk	Black
PTS641965BLXL	Venosan 5001 (Cotton) - Knee - 18-21mmHg - Unisex - C/Toe - XL - Blk	Black
PTS641997	Venosan - Microfibre - Knee - 15-20mmHg - Female - C/Toe - Sml - Bge	Beige
PTS641965MAL	Venosan 5001 (Cotton) - Knee - 18-21mmHg - Unisex - C/Toe - Lge - Marokko	Marokko
PTS641997BEL	Venosan - Microfibre - Knee - 15-20mmHg - Female - C/Toe - Lge - Bge	Beige
PTS641965MAM	Venosan 5001 (Cotton) - Knee - 18-21mmHg - Unisex - C/Toe - Med - Marokko	Marokko
PTS641997BEM	Venosan - Microfibre - Knee - 15-20mmHg - Female - C/Toe - Med - Bge	Beige
PTS641965MAS	Venosan 5001 (Cotton) - Knee - 18-21mmHg - Unisex - C/Toe - Sml - Marokko	Marokko
PTS641997BEXL	Venosan - Microfibre - Knee - 15-20mmHg - Female - C/Toe - XL - Bge	Beige
PTS641965MAXL	Venosan 5001 (Cotton) - Knee - 18-21mmHg - Unisex - C/Toe - XL - Marokko	Marokko
PTS641997BLL	Venosan - Microfibre - Knee - 15-20mmHg - Female - C/Toe - Lge - Blk	Black
PTS641965OBLL	Venosan 5001 (Cotton) - Knee - 18-21mmHg - Unisex - O/Toe - Lge - Blk	Black
PTS641997BLM	Venosan - Microfibre - Knee - 15-20mmHg - Female - C/Toe - Med - Blk	Black
PTS641965OBLM	Venosan 5001 (Cotton) - Knee - 18-21mmHg - Unisex - O/Toe - Med - Blk	Black
PTS641997BLS	Venosan - Microfibre - Knee - 15-20mmHg - Female - C/Toe - Sml - Blk	Black
PTS641965OBLS	Venosan 5001 (Cotton) - Knee - 18-21mmHg - Unisex - O/Toe - Sml - Blk	Black
PTS641997BLXL	Venosan - Microfibre - Knee - 15-20mmHg - Female - C/Toe - XL - Blk	Black

Code	Display Name	Colour
PTS641965OBXL	Venosan 5001 (Cotton) - Knee - 18-21mmHg - Unisex - O/Toe - XL - Blk	Black
PTS641998	Venosan - Microfibre - Knee - 20-30mmHg - Male - C/Toe - Sml - Blk	Black
PTS641965OMAL	Venosan 5001 (Cotton) - Knee - 18-21mmHg - Unisex - O/Toe - Lge - Marokko	Marokko
PTS641998L	Venosan - Microfibre - Knee - 20-30mmHg - Male - C/Toe - Lge - Blk	Black
PTS641965OMAM	Venosan 5001 (Cotton) - Knee - 18-21mmHg - Unisex - O/Toe - Med - Marokko	Marokko
PTS641998M	Venosan - Microfibre - Knee - 20-30mmHg - Male - C/Toe - Med - Blk	Black
PTS641965OMAXL	Venosan 5001 (Cotton) - Knee - 18-21mmHg - Unisex - O/Toe - XL - Marokko	Marokko
PTS641998XL	Venosan - Microfibre - Knee - 20-30mmHg - Male - C/Toe - XL - Blk	Black
PTS641966	Venosan 5002 (Cotton) - Knee - 23-32mmHg - Unisex - O/Toe - Sml - Marokko	Marokko
PTS641966BLL	Venosan 5002 (Cotton) - Knee - 23-32mmHg - Unisex - C/Toe - Lge - Blk	Black
PTS641966BLM	Venosan 5002 (Cotton) - Knee - 23-32mmHg - Unisex - C/Toe - Med - Blk	Black
PTS641966BLS	Venosan 5002 (Cotton) - Knee - 23-32mmHg - Unisex - C/Toe - Sml - Blk	Black
PTS641966BLXL	Venosan 5002 (Cotton) - Knee - 23-32mmHg - Unisex - C/Toe - XL - Blk	Black
PTS641966MAL	Venosan 5002 (Cotton) - Knee - 23-32mmHg - Unisex - C/Toe - Lge - Marokko	Marokko
PTS641966MAM	Venosan 5002 (Cotton) - Knee - 23-32mmHg - Unisex - C/Toe - Med - Marokko	Marokko
PTS641966MAS	Venosan 5002 (Cotton) - Knee - 23-32mmHg - Unisex - C/Toe - Sml - Marokko	Marokko
PTS641966MAXL	Venosan 5002 (Cotton) - Knee - 23-32mmHg - Unisex - C/Toe - XL - Marokko	Marokko
PTS641966OBLL	Venosan 5002 (Cotton) - Knee - 23-32mmHg - Unisex - O/Toe - Lge - Blk	Black

Code	Display Name	Colour
PTS641966OBLM	Venosan 5002 (Cotton) - Knee - 23-32mmHg - Unisex - O/Toe - Med - Blk	Black
PTS641966OBLS	Venosan 5002 (Cotton) - Knee - 23-32mmHg - Unisex - O/Toe - Sml - Blk	Black
PTS641966OBXL	Venosan 5002 (Cotton) - Knee - 23-32mmHg - Unisex - O/Toe - XL - Blk	Black
PTS641966OMAL	Venosan 5002 (Cotton) - Knee - 23-32mmHg - Unisex - O/Toe - Lge - Marokko	Marokko
PTS641966OMAM	Venosan 5002 (Cotton) - Knee - 23-32mmHg - Unisex - O/Toe - Med - Marokko	Marokko
PTS641966OMAXL	Venosan 5002 (Cotton) - Knee - 23-32mmHg - Unisex - O/Toe - XL - Marokko	Marokko