

## Diabpro Sheepskin Hi Temp DiabPro Shoe



- The DiabPro has a narrower sole width compared to other Shear Comfort footwear and has a posted support at the heel. This boot has a slimline design that makes it perfect for everyday use.
- The narrower sole and reduced bulk give the boot a sleek look and are designed for active users who use the boots when walking around the house or outside. Two Velcro straps fasten the boot securely across the top and can be adjusted to accommodate swelling.
- Sheep's wool accommodates all foot shapes by moulding around the foot, adding stability to your everyday walking and protecting delicate skin from friction.

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Clothes and shoes

**Category:** Low risk

Click [here](#) for more information

Code	Display Name	Shoe Size
PTA521910	DiabPro Sheepskin Shoe - Closed-Toe - Extra Extra Small	2 - 3.5
PTA521911	DiabPro Sheepskin Shoe - Closed-Toe - Extra Small	04-May
PTA521912	DiabPro Sheepskin Shoe - Closed-Toe - Small	5.5 - 7
PTA521913	DiabPro Sheepskin Shoe - Closed-Toe - Medium	7.5 - 8.5
PTA521914	DiabPro Sheepskin Shoe - Closed-Toe - Large	09-Oct
PTA521915	DiabPro Sheepskin Shoe - Closed-Toe - Extra Large	10.5 - 12