

Cura - Personal Emergency Transmitter With Lanyard (PET)



- Large rubber button for easy use by elderly patients
- Standard Lanyard for patient to wear around neck
- Works with single use hospital wristbands
- Changeable battery with projected 1-year life
- Splash resistant
- Operating frequency: 433.92MHz
- Operating range: 30 m line of sight
- Power supply: 12V 27A Alkaline Battery
- Battery life: >1 year (depends on usage)

 Contract approved

Support @ Home

Function: Signalling, alarming and localisation

Category: Prescribed

Click [here](#) for more information

Compliance

No applicable compliance standard

cura¹

Code	Display Name	Material
BEA007469	Personal Emergency Transmitter With Lanyard	PET