

OPPO 3187 ARM SLING / IMMOBILISER



- Helps immobilize and support the arm
- Used to comfortably support a sprained or broken arm or for post-surgical support
- Extra strap offers more stabilization to the arm and body Application:
- Brachial nerve injury
- Upper arm fracture
- Shoulder subluxation/dislocation
- Post surgical protection

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name
FBS308520	OPPO 3187 ARM SLING / IMMOBILISER - SML
FBS308510	OPPO 3187 ARM SLING / IMMOBILISER - MED
FBS308500	OPPO 3187 ARM SLING / IMMOBILISER - LGE
FBS308530	OPPO 3187 ARM SLING / IMMOBILISER - XL