

## GelBodies Shin Protector



- Combined with pressure relieving gel, which is found in high-end wheelchair cushions, to offer additional padding to points on the body where pressure sores tend to form
- Improvement of circulation provides even distribution of the body weight and gives patient added overall comfort helping them to relax
- Breathable, stretchy cool-band material
- Sizing: Measure around the calf. Choose the larger size if measurement is between sizes

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Manage tissue integrity

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Calf/Arm Circumference
PTA533326	GelBodies - Shin Protector - Med	25-30cm
PTA533326L	GelBodies - Shin Protector - Lge	33-41cm