

## Rolyan Carve-It Hand Positioning Splint



- Helps anatomically position the hand with mild to moderate finger spasticity
- Design allows use with left or right hand
- Easy to modify by cutting foam with an electric knife for a custom fit
- Foam shaped to fit the hand and comfortable for the patient
- Unlike traditional hand cones, this splint maintains the arches of the hand in an appropriate position
- Flat bottom provides easy positioning in bed or on lap tray, arm tray
- Removable cover can be hand washed and air dried
- For use with strokes, Alzheimer's disease, dementia and other neurological impairments
- For size measure width of 2nd and 5th metacarpophalangeal (MCP) joints and see total width options

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Upper limb orthoses

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name	Total Width
FBS283075	Rolyan Carve-It Hand Positioning Splint - S/M	7.6 cm