

DONJOY PROCARE 79-811 SUPER LITE AFO



- Lightweight and low profile.
- Will fit inside most lace-up shoes.
- Durable, orthopaedic grade thermoplastic may be custom trimmed and reshaped. Application:
- Indicated for foot drop

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name
FBS305240	DONJOY PROCARE 79-811 SUPER LITE AFO - MENS - RIGHT - UNIVERSAL
FBS305250	DONJOY PROCARE 79-811 SUPER LITE AFO - MENS - LEFT - UNIVERSAL
FBS305260	DONJOY PROCARE 79-811 SUPER LITE AFO - WOMENS - RIGHT - UNIVERSAL
FBS305270	DONJOY PROCARE 79-811 SUPER LITE AFO - WOMENS - LEFT - UNIVERSAL