

Allegro - Salsa Standing Transfer Hoist



The SALSA mobile standing hoist facilitates easy and secure transfers from a seated to a supported standing position, providing essential assistance for individuals with limited mobility.



Key Features:

- Electrically operated for effortless use
- Ergonomically designed for patient and caregiver comfort
- Durable construction for long-lasting performance
- Provides a secure and stable standing position

Manufactured by Allegro Concepts, a trusted name in patient lifting equipment, the SALSA Mobile Standing Hoist is part of a comprehensive range that includes various lifters, slings, and accessories, all designed and tested in-house to comply with the industry-standard AS ISO 10535 – 2002.

Warranty Information

- Patient lifter including electronics 2 years
- Accessories on lifter 2 years
- Patient slings 1 year
- Weighing devices 1 year
- Batteries 3 months

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Lift Range	Electric Leg
LSS377410	Allegro - Salsa Standing Transfer Hoist	1010 mm	665 - 1030mm	1190 mm	200 kg	51kg	825-1670mm	True

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Lift Range	Electric Leg
------	--------------	--------------	-------------	--------------	-------------------	-------------	------------	--------------