

## Oppo 4280 Baseball Finger Splint

- Helps stabilise the phalangeal joints in position.
- Malleable aluminium splint can be reshaped to hold finger at an exact degree of flexion or extension.
- Provides firm support while the injury heals.

 Contract approved

### Support @ Home

**Function:** Upper limb orthoses

**Category:** Prescribed

Click [here](#) for more information



Code	Display Name	Total Length
FBS301466	Oppo 4280 Baseball Finger Splint - Small	80 mm
FBS3014662	Oppo 4280 Baseball Finger Splint - Medium	100 mm