

Sheepskin Hi Temp Classic Mattress Overlay



- A full covered medical sheepskin bed overlay, designed to cover head to toe
- Suitable for long term patients where shear and friction are a concern

Warranty Information

- 1 Year

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Not on S@H

Category: Not on S@H

Click [here](#) for more information

NDIS

Item Code:

05_043306004_0103_1_2

Power Bed – Composite:

05_423_0103_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Total Length	Total Width
PTA521400	Sheepskin Mattress Overlay	1250 mm	750 mm