

ROHO - LTV



permobil

ROHO

The LTV Seat Cushion is a non-medical cushion application for anyone who sits for extended periods of time. The LTV Seat Cushion is lightweight and easy to fold, making this the perfect cushion to take anywhere. Shape Fitting Technology® provides unsurpassed comfort. The LTV is made of interconnected air cells that allow air to slowly transfer from chamber to chamber and evenly distribute body weight pressure.

- Appropriate for everyday comfort and pressure injury prevention
- Interconnected air-filled cells to redistribute body weight and minimise pressure
- Padded cover or optional synthetic leather cover adds an extra level of comfort
- Lightweight design and air valve make it easy to fold and great for travel
- Height adjustable 1 - 2"

Warranty Information

- 1 Year

✓ Contract approved

Support @ Home

Function: Manage tissue integrity

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_181006190_0105_1_2

Pressure Care Cushion –

Composite:05_429_0103_1_2

Click [here](#) for more information

Compliance

EN-ISO 10993-10 (1999-11) / EN 1021-1 and EN 1021-2

Code	Display Name
PTC585200	ROHO LTV Cushion