

Action 4NG



- Invacare Action 4NG is comfortable, foldable, and lightweight, making it suited for everyday use.
- Cross-compatibility allows users and advisors to update and refurbish their wheelchairs with other models in the Action NG family.
- Flip back armrests.
- Step-less adjustable backrest angle.
- Swing away/removable footplates.

 Contract approved

Support @ Home

Function: Manual wheelchairs

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height
MWC427050	Wheelchair - Action 4NG (As per script)	160 kg	14 kg	450 - 500 mm	505 - 605 mm	510 mm