

Teal Tolero High Back Chair



Tolero is a metal patient chair with full feature set, including removable seat, transport wheels, height adjustable legs, slim contoured profile and infection management features. Designed for use in hospital and healthcare environments.



- High back height adjustable day chair with 100mm seat height variation
- Slim, contoured ergonomically designed posture back for head/lumbar support
- Strong, durable high-grade 40mm tubular metal frame creates a modern aesthetic
- Anti-bacterial easy-clean coated finish for infection control management
- Removable, depth adjustable seat as standard with sacral gap for cleaning/access
- Integrated transport wheels
- 100% recyclable material
- Tested and compliant to BS EN 16139:2013

Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height
CHP231300	Teal Tolero High Back Chair - (Choose colour)	720 mm	608 mm	1195 - 1295 mm	160 kg	37 kg	445 mm	470 mm	460 - 560 mm