

Teal Tolero High Back Chair



Tolero is a metal patient chair with full feature set, including removable seat, transport wheels, height adjustable legs, slim contoured profile and infection management features. Designed for use in hospital and healthcare environments.

- High back height adjustable day chair with 100mm seat height variation
- Slim, contoured ergonomically designed posture back for head/lumbar support
- Strong, durable high-grade 40mm tubular metal frame creates a modern aesthetic
- Anti-bacterial easy-clean coated finish for infection control management
- Removable, depth adjustable seat as standard with sacral gap for cleaning/access
- Integrated transport wheels
- 100% recyclable material
- Tested and compliant to BS EN 16139:2013

Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Unit Weight	Seat Depth	Seat Width	Seat To Floor Height	Seat To Footrest
CHP231300	Teal Tolero High Back Chair - (Choose colour)	720 mm	608 mm	1195 - 1295 mm	190 kg	37 kg	445 mm	470 mm	460 mm	460 mm