

Theraband Tubing



- TheraBand Exercise Tubing is a natural rubber tubing available in 6 colour-coded levels of resistance
- Improves strength, range of motion and cooperation of muscle groups

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Not on S@H

Category: Not on S@H

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Colour
FPP310800	7.5m Special Heavy TheraBand Tubing	Black
FPP309450	30m Super TheraBand Tubing	Silver
FPP309500	30m Special Heavy TheraBand Tubing	Black
FPP309600	30m Extra Heavy TheraBand Tubing	Blue
FPP309700	30m Heavy TheraBand Tubing	Green
FPP309800	30m MediumTheraBand Tubing	Red
FPP310850	7.5m Super Heavy TheraBand Tubing	Silver

Code	Display Name	Colour
FPP310700	7.5m Extra Heavy TheraBand Tubing	Blue
FPP310600	7.5m Heavy Theraband Tubing	Green
FPP310500	7.5m Medium Theraband tubing	Red
FPP310400	7.5m Low TheraBand Tubing	Yellow
FPP310000	Theraband - Exercise Handles	-