

Multistrap



- The MultiStrap lift aids are lifting accessories that complete the menu of lifting and repositioning aids in Liko's wide product range
- The MultiStrap aid was developed to meet a variety of needs in care and can be used to hold arms, legs and other parts of the body during treatment, washing or other routine caregiving activities
- Together with a patient lift, the MultiStrap lift aids ease many heavy lifting tasks such as wound care, hygiene procedures, catheterisation, examinations, placement of bedpans, dressing/undressing and more
- MultiStrap lift aids may also be used to turn a patient in bed
- They are available in Standard 10 cm, Wide 25 cm and Ultra 40 cm models, and can be used both individually and in pairs
- Also available in a Solo version for single-patient use to facilitate infection control compliance
- FixStrap accessory is available to hold the MultiStrap in place during turning



Compliance

AS/NZS ISO 10535:2011

Warranty Information

- 1 year

Code	Display Name	Safe Working Load
LSS395300	MultiStrap Std 10cm - Liko - Poly	200kg
LSS395310	MultiStrap Wide 25cm - Liko - Poly	200kg
LSS395320	MultiStrap Ultra - Liko - Poly	200kg

