

## LIKO ORIGINAL HIGHBACK SLING (MOD. 200/210)



- This sling offers a reclined seating position for patients with poor head and trunk stability
- Separated leg supports ensure that the patient does not sit on the sling after use
- Easily adaptable to the patient with comfortable full body support
- Can be used for floor transfers
- Available in Polyester or Mesh

### Warranty Information

- 1 Year



 Contract approved

**ndis**

**Item Code:**

05\_123603164\_0105\_1\_2

**Hoist –**

**Composite:**05\_420\_0105\_1\_2

Click [here](#) for more information

**Compliance**

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS394650	Original Highback Sling - Reinforced Leg & Belt - Polyester - XLARGE	500 kg	-
LSS394660	Original Highback Sling - Reinforced Leg & Belt - Mesh - XSMALL	200 kg	Mesh (Loops)
LSS394670	Original Highback Sling - Reinforced Leg & Belt - Mesh - SMALL	200kg	Mesh (Loops)
LSS394680	Original Highback Sling - Reinforced Leg & Belt - Mesh - MEDIUM	200 kg	Mesh (Loops)
LSS394690	Original Highback Sling - Reinforced Leg & Belt - Mesh - LARGE	200 kg	Mesh (Loops)
LSS394700	Original Highback Sling - Reinforced Leg & Belt - Mesh - XLARGE	500 kg	Mesh (Loops)
LSS394610	Original Highback Sling - Reinforced Leg & Belt - Polyester - XSMALL	200 kg	Polyester (Loops)

Code	Display Name	Safe Working Load	Material
LSS394620	Original Highback Sling - Reinforced Leg & Belt - Polyester - SMALL	200 kg	Polyester (Loops)
LSS394630	Original Highback Sling - Reinforced Leg & Belt - Polyester - MEDIUM	200 kg	Polyester (Loops)
LSS394640	Original Highback Sling - Reinforced Leg & Belt - Polyester - LARGE	200 kg	Polyester (Loops)
LSS394675	Original Highback Sling - Reinforced Leg Support - Net Polyester - MEDIUM SLIM	200 kg	Net Polyester
LSS394685	Original Highback Sling - Reinforced Leg Support - Polyester - MEDIUM SLIM	200 kg	Polyester