

Aspire Bath Transfer Bench - Padded



- Padded seat & backrest provides a cushioned surface for comfort
- Highly resilient vinyl covering
- Fully height adjustable with a side arm for secure transfers
- Adjustable backrest to suit right or left facing tubs
- Suction cup foot pieces secure the bench to the tub floor
- Lightweight aluminium frame construction
- **Aspire Extension Legs - [Set of 4](#)**

Each set includes:

- 2 x suction tips (for use on smooth surfaces)
- 2 x standard tips (for general use)

✓ Contracted DVA product

✓ Contract approved

ndis

Item Code:

05_123121160_0103_1_2

Shower Support –

Composite:05_415_0103_1_2

Click [here](#) for more information

Warranty Information

- 1 Year



Compliance

AS/NZS 3973:2009

Code	Display Name	RRP	Box Quantity	Seat Depth	Seat Width	Seat To Footrest	Safe Working Load	Unit Weight
BTS110900	Aspire Bath Transfer Bench - Padded	292	1	405 mm	600 mm	455 - 555 mm	181 kg	4.9 kg
BTS108322	Extension Legs - 100mm - w/2 Suction & 2 Std Tips - Suits BTS110900	60	1	-	-	-	-	-