

## Aspire Bath Transfer Bench - Padded



- Padded seat & backrest provides a cushioned surface for comfort
- Highly resilient vinyl covering
- Fully height adjustable with a side arm for secure transfers
- Adjustable backrest to suit right or left facing tubs
- Suction cup foot pieces secure the bench to the tub floor
- Lightweight aluminium frame construction
- **Aspire Extension Legs** - [Set of 4](#)

Each set includes:

- 2 x suction tips (for use on smooth surfaces)
- 2 x standard tips (for general use)

✓ Contracted DVA product

✓ Contract approved

**ndis**

**Item Code:**

05\_123121160\_0103\_1\_2

**Shower Support –**

**Composite:**05\_415\_0103\_1\_2

Click [here](#) for more information

### Warranty Information

- 1 Year



**Compliance**

AS/NZS 3973:2009

Code	Display Name	Seat Depth	Seat Width	Seat To Footrest	Safe Working Load	Unit Weight
BTS110900	Aspire Bath Transfer Bench - Padded	405 mm	600 mm	455 - 555 mm	181 kg	4.9 kg
BTS108322	Extension Legs - 100mm - w/2 Suction & 2 Std Tips - Suits BTS110900	-	-	-	-	-