

Freedom Height Adjustable Bed Grab Bar - JBA506



- Designed to assist the patient into a seated position
- Padded hand grip with curves for safety and comfort
- Height adjustable
- Can be used on either side of the bed
- Powder-coated aluminium frame
- Fits securely under mattress

✓ Contracted DVA product

✓ Contract approved

Compliance

AS EN 12182:2015

Warranty Information

- 1 Year

freedom
HEALTHCARE

Code	Display Name	RRP	Box Quantity	Total Length	Total Width	Total Height	Material	Unit Weight
BEA831100	Freedom Height Adjustable Bed Grab Bar	49.5	1	550 mm	830 mm	420-550 mm	Aluminium	1.9 kg