

## Freedom Height Adjustable Bed Grab Bar - JBA506



- Designed to assist the patient into a seated position
- Padded hand grip with curves for safety and comfort
- Height adjustable
- Can be used on either side of the bed
- Powder-coated aluminium frame
- Fits securely under mattress

### Warranty Information

- 1 Year

**freedom**  
HEALTHCARE

✓ Contracted DVA product

✓ Contract approved

### Support @ Home

**Function:** Beds and bed equipment

**Category:** Prescribed

Click [here](#) for more information

### Compliance

AS EN 12182:2015

Code	Display Name	Total Length	Total Width	Total Height	Material	Unit Weight
BEA831100	Freedom Height Adjustable Bed Grab Bar	550 mm	830 mm	420-550 mm	Aluminium	1.9 kg