

Freedom Height Adjustable Bed Grab Bar - JBA506



- Designed to assist the patient into a seated position
- Padded hand grip with curves for safety and comfort
- Height adjustable
- · Can be used on either side of the bed
- Powder-coated aluminium frame
- Fits securely under mattress

Warranty Information

• 1 Year



Code	Display Name	Total Length	Total Width	Total Height	Material	Unit Weight
BEA831100 Fre	eedom Height Adjustable Bed Grab Bar	550 mm	830 mm	420-550 mm	Aluminium	1.9 kg

Ocontracted DVA product

⊘Contract approved

Compliance

AS EN 12182:2015