

## Invacare General Purpose Sling with Head Support



- This sling fits easily to a patient in bed or a chair, providing excellent upper body and head support
- Recommended for lifting from the floor and general lifting of the less secure patient



**ndis**

**Item Code:**

05\_123603164\_0105\_1\_2

**Sling – Composite:**

05\_421\_0105\_1\_2

Click [here](#) for more information

### Warranty Information

- 1 Year



### Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS374300	Invacare General Purpose Poly with Head Support - Extra Small	205 kg	Polyester
LSS373200	Invacare General Purpose Poly with Head Support - Medium	205 kg	Polyester
LSS373300	Invacare General Purpose Poly with Head Support - Large	300 kg	Polyester
LSS374310	Invacare General Purpose Poly with Head Support - Extra Large	300 kg	Polyester
LSS373700	Invacare General Purpose Mesh with Head Support - Small	205 kg	Mesh
LSS373704	Invacare General Purpose Mesh with Head Support - Medium	205 kg	Mesh
LSS373706	Invacare General Purpose Mesh with Head Support - Large	300 kg	Mesh

Code	Display Name	Safe Working Load	Material
LSS374910	Invacare General Purpose Mesh with Head Support - Extra Large	300 kg	Mesh