

Invacare General Purpose Sling with Head Support



- This sling fits easily to a patient in bed or a chair, providing excellent upper body and head support
- Recommended for lifting from the floor and general lifting of the less secure patient

Warranty Information

- 1 Year



ndis

Item Code:

05_123603164_0105_1_2

Sling – Composite:

05_421_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

| Code | Display Name | Safe Working Load | Material |
|-----------|---|-------------------|-----------|
| LSS374300 | Invacare General Purpose Poly with Head Support - Extra Small | 205 kg | Polyester |
| LSS374909 | Invacare General Purpose Poly with Head Support - Small | 205 kg | Polyester |
| LSS373200 | Invacare General Purpose Poly with Head Support - Medium | 205 kg | Polyester |
| LSS373300 | Invacare General Purpose Poly with Head Support - Large | 300 kg | Polyester |
| LSS374310 | Invacare General Purpose Poly with Head Support - Extra Large | 300 kg | Polyester |
| LSS373700 | Invacare General Purpose Mesh with Head Support - Small | 205 kg | Mesh |
| LSS373704 | Invacare General Purpose Mesh with Head Support - Medium | 205 kg | Mesh |

| Code | Display Name | Safe Working Load | Material |
|-----------|---|-------------------|----------|
| LSS373706 | Invacare General Purpose Mesh with Head Support - Large | 300 kg | Mesh |
| LSS374910 | Invacare General Purpose Mesh with Head Support - Extra Large | 300 kg | Mesh |