

MoliCare Premium Form



- Pads for moderate to severe incontinence
- Shaped to fit perfectly to the body
- 3 layer absorbent core keeps skin dry
- Maintains the skin pH and reduces the risk of incontinence related skin irritation
- Inner leakage cuffs for maximum leakage protection
- Odour neutralisers for high level confidence
- · Soft textile like materials for ultimate wearing comfort and discretion
- Recommended to be worn with MoliCare Premium Fixpants



Compliance

No applicable compliance standard

Code	Display Name	Capacity
------	--------------	----------

CON252600 MoliCare Form - Normal Plus (Pkt 30) 4 Drops (1483ml)