

## Panamatic Mini



- The Panamatic Mini - Flusher Thermal Sanitiser is the cold water cousin of the Optima 2 and is designed for low-medium workload situations
- Improved Wash Performance (Cold water) - An increased number of more efficient cleaning nozzles and optimum water jet pressure result in performance levels that make the Mini outstanding for its size
- Reduced Cycle-Time - Better washing and temperature control has enabled improved overall cycle-time for quicker turnaround
- More Efficient Water and Power - With improved performance and reduced cycle time, both power and water consumption have been reduced. Cold water flush / wash - steam cycle - cool-down rinse
- Hands-free Opening and Cycle Start - A foot operated opening mechanism and hands-free close and cycle start means direct contact with the machine has been reduced to a minimum, helping to reduce potential cross infection
- MicrobeSafe Antimicrobial Lid - DDC™ unique lid incorporates Silver Nanotechnology to help stop bacterial growth and reduce cross infection
- Large Chamber Capacity and Easier Loading - The wash chamber and larger opening allow easier loading of a wide range of utensils. The chamber can process 2 bedpans and 2 urine bottles per load or various combinations
- Electrical supply: Single Phase 230 V, 13 AMP
- Electrical Output / Input: Pump - 0.425 kW
- Water supply: 15 mm Cold water 0.5 – 10 bar
- Soil connection: 75 mm connections supplied
- Operation cycle time: Approx 6 minutes / 30 ml Ecowash
- Steam Disinfection: Exceeds Ao 700, 85 degrees x 75 secs



### Compliance

ISO 15883/AS 4187:2019

## Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Capacity
UPM100290	DDC Bedpan Washer Sanitiser - Mini - 1 Phase	500 mm	600 mm	1 bedpan & 2 urine bottles / cycle