

## Aspire Viking Forearm Walker



The Viking Forearm Walker is designed for hospital and clinical use, providing stable, comfortable, and fully adjustable support for upright walking. Its strong aluminium frame, padded arm supports, and easy-to-use hydraulic height adjustment make it ideal for rehabilitation and assisted mobility in care settings.

- **Hydraulic Height Adjustment:** Lever-operated system makes it easy to raise or lower the arm supports smoothly and safely, helping users of different heights
- **Adjustable Padded Arm Supports:** Soft armrests offer comfortable forearm support, for users with limited grip strength. The width can be adjusted to suit individual body shapes and provide added stability
- **Width-adjustable Base:** Improves manoeuvrability in tight spaces and around clinical equipment
- **Angle-adjustable Soft Handgrips:** Supports a natural wrist position, making it easier and more comfortable for users to maintain control
- **Oxygen Tank Holder:** Included as standard, which can be fitted or removed to suit the users needs
- **Optional Accessories:** Available to enhance functionality

✓ Contracted DVA product

✓ Contract approved

### Support @ Home

**Function:** Walking, manipulated by both arms

**Category:** Prescribed

Click [here](#) for more information

### Compliance

ISO 11199-2:2005

## Warranty Information

- 1 Year

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Colour
WAF699130	Aspire Viking Forearm Walker	850 mm	730 - 890 mm	900 - 1325 mm	200 kg	Aluminium	Blue

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Colour
GSP650520	Spare Part - Rear Lockable Castor Only -	-	-	-	-	-	-
WAA648340	Bag - Viking Forearm Walker	-	-	-	-	-	-
WAA648350	Bag Holder - Viking Forearm Walker	-	-	-	-	-	-
WCA871190	Iv pole - Viking Forearm Walker	-	-	-	-	-	-