

Body Assist N71 Thermal Shoulder Support



- Provides soft tissue support & compression
- Increases elasticity of muscles, ligaments & tendons
- Assists in prevention of injury or re-injury
- Elastic support strap fits comfortably under the arm pit & assists to generate better than average compression & support to the shoulder
- Can be worn while playing sports
- Multi-directional compression assists to reduce soft tissue swelling
- Promotes greater blood circulation, provides a more speedy recovery
- Contoured panels fit most body shapes
- Completely washable & hygienic
- Stabiliser strap is optional
- Note: Measure circumference around chest, around outside of arm on affected side

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Compliance

No applicable compliance standard

| Code | Display Name | Chest Circumference |
|-----------|--|---------------------|
| FBS281398 | Body Assist N71 Thermal Shoulder Support - Left - X Sml | 80 - 90 cm |
| FBS281399 | Body Assist N71 Thermal Shoulder Support - Left - Sml | 90 - 100 cm |
| FBS281400 | Body Assist N71 Thermal Shoulder Support - Left - Med | 100 - 110 cm |
| FBS281401 | Body Assist N71 Thermal Shoulder Support - Left - Lge | 110 - 120 cm |
| FBS281402 | Body Assist N71 Thermal Shoulder Support - Left - XL | 120 - 130 cm |
| FBS281403 | Body Assist N71 Thermal Shoulder Support - Right - X Sml | 80 - 90 cm |

| Code | Display Name | Chest Circumference |
|-------------|--|----------------------------|
| FBS281405 | Body Assist N71 Thermal Shoulder Support - Right - Sml | 90 - 100 cm |
| FBS281406 | Body Assist N71 Thermal Shoulder Support - Right - Med | 100 - 110 cm |
| FBS281407 | Body Assist N71 Thermal Shoulder Support - Right - Lge | 110 - 120 cm |
| FBS281408 | Body Assist N71 Thermal Shoulder Support - Right - XL | 120 - 130 cm |