

Body Assist N71 Thermal Shoulder Support



- Provides soft tissue support & compression
- Increases elasticity of muscles, ligaments & tendons
- Assists in prevention of injury or re-injury
- Elastic support strap fits comfortably under the arm pit & assists to generate better than average compression & support to the shoulder
- Can be worn while playing sports
- Multi-directional compression assists to reduce soft tissue swelling
- Promotes greater blood circulation, provides a more speedy recovery
- Contoured panels fit most body shapes
- Completely washable & hygienic
- Stabiliser strap is optional
- Note: Measure circumference around chest, around outside of arm on affected side

✓ Contracted DVA product

✓ Contract approved

Compliance

No applicable compliance standard

Code	Display Name	Chest Circumference
FBS281398	Body Assist N71 Thermal Shoulder Support - Left - X Sml	80 - 90 cm
FBS281399	Body Assist N71 Thermal Shoulder Support - Left - Sml	90 - 100 cm
FBS281400	Body Assist N71 Thermal Shoulder Support - Left - Med	100 - 110 cm
FBS281401	Body Assist N71 Thermal Shoulder Support - Left - Lge	110 - 120 cm
FBS281402	Body Assist N71 Thermal Shoulder Support - Left - XL	120 - 130 cm
FBS281403	Body Assist N71 Thermal Shoulder Support - Right - X Sml	80 - 90 cm

Code	Display Name	Chest Circumference
FBS281405	Body Assist N71 Thermal Shoulder Support - Right - Sml	90 - 100 cm
FBS281406	Body Assist N71 Thermal Shoulder Support - Right - Med	100 - 110 cm
FBS281407	Body Assist N71 Thermal Shoulder Support - Right - Lge	110 - 120 cm
FBS281408	Body Assist N71 Thermal Shoulder Support - Right - XL	120 - 130 cm