

Oppo 1083 Wrist Wrap with Thumb Loop



- Provides comfortable compressions and support for weak or stressed wrists
- Helps treatment and prevention of injury
- Retains warmth, which may aid and speed up recovery period
- Universal size

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Upper limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Wrist Circumference
FBS301368	Oppo - Wrist Wrap - Thumb Loop - Neoprene - (UNIVERSAL)	120 - 230 mm