

JAY - J2 & J2 Tall Backrest




J2 Tall

- Designed for the client who has minimal to moderate postural needs that may change over time
- The lightweight contoured aluminium shell and contoured foam insert along with 15° of tilt and height adjustability replace the standard wheelchair back upholstery
- J2 provides excellent posterior pelvis stabilisation and posterior thoracic support
- The J2 also has multiple compatible options such as swing-away lateral supports, offset hardware kits and fluid pads to reduce spinal pressure points
- Contoured foam insert distributes pressure over wider area to increase comfort and reduce the risk of skin breakdown

Warranty Information

- 2 Years Backrest
- 6 Months Cover

 Contract approved

ndis

Item Code:

05_181003192_0105_1_2

Wheelchair Manual – Composite:

Wheelchair Powered –

Composite:05_417_0105_1_2

05_419_0105_1_2

Click [here](#) for more information

Compliance

No applicable compliance standard

Code	Display Name	Seat Width	Backrest Width	Backrest Height	Safe Working Load
PTB532260	JAY J2 - 12in	11 - 13 in	12 in	12.5 in	113 kg
PTB532265	JAY J2 - 13in	12 - 14 in	13 in	13 in	113 kg
PTB532270	JAY J2 - 14in	13 - 15 in	14 in	14.5 in	113 kg
PTB532280	JAY J2 Tall - 14in	13 - 15 in	14 in	21 in	113 kg
PTB532300	JAY J2 Tall - 15in	14 - 16 in	15 in	21 in	113 kg
PTB532310	JAY J2 - 16in	15 - 17 in	16 in	16.5 in	113 kg
PTB532320	JAY J2 Tall - 16in	15 - 17 in	16 in	22 in	113 kg
PTB532330	JAY J2 - 17in	16 - 18 in	17 in	16.5 in	113 kg
PTB532340	JAY J2 Tall - 17in	16 - 18 in	17 in	22 in	113 kg
PTB532350	JAY J2 - 18in	17 - 19 in	18 in	16.5 in	113 kg
PTB532360	JAY J2 Tall - 18in	17 - 19 in	18 in	22 in	113 kg
PTB532370	JAY J2 - 19in	18 - 20 in	19 in	16.5 in	113 kg
PTB532380	JAY J2 Tall - 19in	18 - 20 in	19 in	23 in	113 kg
PTB532390	JAY J2 - 20in	19 - 21 in	20 in	17.5 in	113 kg
PTB532400	JAY J2 Tall - 20in	19 - 21 in	20 in	23 in	113 kg