

## LifeSmart Blood Glucose Monitor - Bluetooth



### Smart Blood Glucose plus B-Ketone Monitoring System

- Unique blood glucose monitor accurately measures blood glucose, hematocrit and ketone levels.
- Wide hematocrit range (0-70%) enables it to be used on infants to old age adults.
- Actively connects to your smartphone (Android or Apple) allowing you to transfer, view and analyze your meaningful data in real time.
- The free app, gives you the flexibility to share your data with your doctor, dietician or family members.
- Ketone warning.
- 7, 14, 21, 28, 60, 90 day average
- 1000 memory sets with Date/Time/AC/PC
- LifeSmart Model Number: LS946N

✓ Contracted DVA product

