

## MEDI R.210 THUMB SUPPORT



- Used for the orthotic fitting of the thumb area and only on intact skin.
- Skin-friendly material
- Mouldable aluminium splint
- For the left and right hand

✔ Contracted DVA product

✔ Contract approved

### Support @ Home

**Function:** Upper limb orthoses

**Category:** Prescribed

Click [here](#) for more information

Code	Display Name
FBS307510	MEDI R.210 THUMB SUPPORT - LEFT - MED
FBS307520	MEDI R.210 THUMB SUPPORT - LEFT - LGE
FBS307530	MEDI R.210 THUMB SUPPORT - LEFT - XL
FBS307540	MEDI R.210 THUMB SUPPORT - RIGHT - SML
FBS307550	MEDI R.210 THUMB SUPPORT - RIGHT - MED
FBS307560	MEDI R.210 THUMB SUPPORT - RIGHT - LGE
FBS307570	MEDI R.210 THUMB SUPPORT - RIGHT - XL

**Code**

**Display Name**

FBS307500 MEDI R.210 THUMB SUPPORT - LEFT - SML