

Body Assist 41S Cross-Cut Elastic Knee Brace



- Anatomically shaped by a cross-cut design that minimises irritation & pressure behind the knee
- Provides orthopaedic support for minor knee joint strains & sprains
- Removable medical support spring stays on sides of the knee joint to provide medio-lateral support
- Brace retains some body warmth & can increase blood circulation
- Lightweight, ventilated elastic for all day comfort

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lower limb orthoses

Category: Prescribed

Click [here](#) for more information

Code	Display Name	Knee Circumference
FBS284180	BodyAssist - Knee Brace - Elastic Cross-Cut w/ Rods - Sml	up to 320 mm
FBS284190	BodyAssist - Knee Brace - Elastic Cross-Cut w/ Rods - Med	320 - 370 mm
FBS284200	BodyAssist - Knee Brace - Elastic Cross-Cut w/ Rods - Lge	370 - 420 mm
FBS284210	BodyAssist - Knee Brace - Elastic Cross-Cut w/ Rods - XL	420 - 470 mm
FBS284215	BodyAssist - Knee Brace - Elastic Cross-Cut w/ Rods - XXL	470 mm+