

## Sheepskin Hi Temp Medical Overlay



- Medical grade sheepskin is a versatile product that can be used in the bed, on a chair or on a car seat
- Assists in pressure reduction, will absorb water, transfer heat away from the body and reduce shear and friction

### Warranty Information

- 1 Year

- ✓ Contracted DVA product
- ✓ Contract approved

### ndis

#### Item Code:

05\_043306004\_0103\_1\_2

#### Power Bed – Composite:

05\_423\_0103\_1\_2

Click [here](#) for more information

### Compliance

No applicable compliance standard

Code	Display Name	RRP	Box Quantity	Total Length	Total Width
PTA521500	Medical Sheepskin Overlay	259	1	1000 mm	650 mm