

Oxford Advance Folding Patient Lifter



- Mobile patient lifter for transfer to a seated or lying position
- Compact design can be folded easily for storage or transportation
- Epoxy powder coated frame
- Single speed handset control
- Emergency lowering
- Over-sized push handles for excellent manoeuvrability
- Manual foot operated leg gating and tapered leg design for close access
- Suits a large range of patient slings
- Supplied with battery and charger
- Item Weight: 31.7 kg

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

Warranty Information

- 2 Years



NDIS

Item Code:

05_123603161_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Safe Working Load	Material	Unit Weight	Lift Range	Leg Gait Range
LSS377300	Advance Folding Patient Lifter	1300 mm	550 mm	1360 - 1860 mm	155 kg	Aluminium	31.7 kg	350 - 1860 mm	570 - 1000 mm