

K Care Shower Commode Swing Away Footrest



- Assists transportation of the user for showering and can be positioned over the toilet for toileting
- Stainless steel frame for corrosion resistance
- Fully height adjustable legrests to suit user needs
- Padded or plastic flip up toilet seat available
- Mounted on 4 castors (2 are locking for stability)
- Swing aways legrests and arms can be moved out of the way to assist transfers
- Note: Pictured with an optional closed padded commode seat (BTC075500)
- **SEATS, PANS & PAN-CARRIERS ARE SOLD AS SEPARATE ITEMS TO THE COMMODE**



Warranty Information

- 1 Year