

Sphygmomanometer Digital Auto



- Fully automatic providing accurate and comfortable upper arm blood pressure measurement with Enhanced IntelliSense Technology.
- Personalised inflation for maximum comfort
- Quick deflation releases valve for speedy measurement
- Stores 30 Measurements in Memory
- Body Movement Detection
- Irregular Heartbeat Detection
- Hypertension Indicator

Warranty Information

- 1 Year

✓ Contracted DVA product

✓ Contract approved

Support @ Home

Function: Not on S@H

Category: Not on S@H

Click [here](#) for more information

Code

Display Name

GMD300330 Sphygmomanometer Digital Auto