

Oxford Pivot Frames



- The 4-point positioning cradle uses the Securi3 safety clip system and adjustment is made by rotating the cradle handle either down or up to obtain an upright or reclined position for the patient
- Both manual and powered cradles are available
- Suits for Oxford Clip Slings
- LSS377395 suits Advance Folding Patient Lifter (LSS377300)
- LSS377396, LSS377397, LSS377385, and LSS377386 suits Presence Patient Lifter (LSS377325) and Stature Patient Lifter (LSS377373)
- LSS377388 suits Presence Patient Lifter (LSS377325)

Warranty Information

- 2 Years



 Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load
LSS377396	Pivot Frame - 4 Point - Stature & Presence Patient Lifters	227 kg
LSS377385	Pivot Frame with Weigh Scale - 4 Point - Oxford Pro Patient Lifters	227 kg