

Venosan - 'Silverline'



- Everyday Comfort - 20 - 30mmHg compression SILVERLINE offers long-lasting protection from foot odour and athlete's foot
- X-Static, the silver fibre, promotes hygiene and anti-microbial properties
- Skin remains deodorised and free of bacteria and fungus, particularly effective in protecting against foot odour and athlete's foot
- Cool and dry in Summer and warm in Winter, SILVERLINE fibres wick moisture away from the skin for comfort in warm weather
- X-Static, the silver fibre, conducts heat in Winter to retain warmth
- The silver in X-Static yarns is an excellent conductor and eliminates static electricity
- Specifically designed for enhanced comfort and fit with knitted top band
- Elimination of bacteria and fungus associated with foot odour
- Excellent regulator of heat

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Code	Display Name	Colour
PTS650335	Venosan - 'Silverline' - Knee - 20-30mmHg - Female - C/Toe - XXLge - Bge	Beige
PTS641994	Venosan - Silveline - Knee - 20-30mmHg - Male - C/Toe - Sml - Blk	Black
PTS641994L	Venosan - Silveline - Knee - 20-30mmHg - Male - C/Toe - Lge - Blk	Black
PTS641994M	Venosan - Silveline - Knee - 20-30mmHg - Male - C/Toe - Med - Blk	Black
PTS641994XL	Venosan - Silveline - Knee - 20-30mmHg - Male - C/Toe - XL - Blk	Black
PTS641995	Venosan - Silveline - Knee - 20-30mmHg - Female - C/Toe - Sml - Bge	Beige

Code	Display Name	Colour
PTS641995BEL	Venosan - Silveline - Knee - 20-30mmHg - Female - C/Toe - Lge - Bge	Beige
PTS641995BEM	Venosan - Silveline - Knee - 20-30mmHg - Female - C/Toe - Med - Bge	Beige
PTS641995BEXL	Venosan - Silveline - Knee - 20-30mmHg - Female - C/Toe - XL - Bge	Beige
PTS641995BLL	Venosan - Silveline - Knee - 20-30mmHg - Female - C/Toe - Lge - Blk	Black
PTS641995BLM	Venosan - Silveline - Knee - 20-30mmHg - Female - C/Toe - Med - Blk	Black
PTS641995BLS	Venosan - Silveline - Knee - 20-30mmHg - Female - C/Toe - Sml - Blk	Black
PTS641995BLXL	Venosan - Silveline - Knee - 20-30mmHg - Female - C/Toe - XL - Blk	Black
FBS310160	Venosan - Silverline - Knee - 20-30Mmhg - Male - C/Toe - XXLGE - Blk	Black
FBS310170	Venosan - Silverline - Knee - 20-30Mmhg - Female - C/Toe - XXLGE - Blk	Black