

Austin Bedside Cabinet



- The Austin Bedside Cabinet delivers simple, understated styling and generous storage
- Raised safety edging on the cabinet top prevents items from sliding off
- Four drawers help to reduce clutter and keep essential items at hand
- Rear towel rail offers a practical and tidy solution for compact living
- Slimline handles and smooth drawer action offer function and practicality
- A versatile bedside cabinet that suits most settings
- Laminate finish is durable and easy to care for
- Select from our range of laminate colour options
- Optional hutch provides added shelf space along with practical storage for books and magazines



| Code | Display Name | Total Length | Total Width | Total Height |
|-----------|---|--------------|-------------|--------------|
| BEA826300 | Austin Bedside Cabinet - Laminate | 450 mm | 465 mm | 810 mm |
| BEA826305 | Austin Bedside Cabinet w/Hutch - Laminate | 450 mm | 465 mm | 1020 mm |