

## NuboMed M3100-F Medication Auxiliary Cabinet



NuboMed

- **User Interaction:** Features a 15.6" touch screen monitor with a Capacitive Fingerprint Reader, optional facial recognition (with mask recognition capability), Staff RFID card reader, Built-in monitor speaker, and a Barcode Scanner (1D/2D barcodes)
- **Cabinet Material:** The cabinet is composed of sheet metal and high-density foam insulation, sheet metal is painted matte white
- **Cabinet Control Lock:** Intelligent Lock Control; Authorised user access only by login in
- **Cabinet Lighting:** Auto lighting on/off when door is opened/closed
- **Cabinet Door:** Double-layer electrically heated tempered glass (safety shatterproof), defrost and anti-fog function
- **Smart Drawer:** Up to 5 layers of Smart Shelves for medication storage. Each steel shelf has 4x 1U1H medicine bins and an intelligent monitoring response time of 0.5s or less
- **Proximity Scanner:**
  - Recognition distance: <50cm
  - Recognition speed: <2 seconds
- **Guiding Light:** Each bin has a RGB indicator light as a guide to direct operation
- **Smart Colour Bin Display:** The mini colour screen displays the name, specification, and instructions if available
- **Temperature and Humidity Monitoring System:** Smart real-time detection of cabinet temperature and humidity along with a digitalised LED panel display on the cabinet
- **Temperature and Humidity Alarm:** Smart detection of temperature and humidity with auto alarm if exceeding the range
- **Ideal For:** Outpatient Pharmacy, Emergency Room, Pharmacy, Central



Pharmacy and Clinical Ward

## Warranty Information

- 2 Years

Code	Display Name	Total Length	Total Width	Total Height	Unit Weight
WAT908564	NuboMed M3100-F Medication Auxiliary Cabinet	752 mm	750 mm	1967 mm	200 kg