

Spex Vigour - Mid Active Back Support



- Lightweight and designed for efficient propelling, with scapular subtractions for unimpeded arm function.
- A solid back support with the ability to self-propel
- Improved spinal and pelvic positioning whilst maintaining propelling efficiency
- Minimal to moderate postural support
- In combination with lateral supports offers exceptional trunk and thoracic stability
- For clients that have mild torso asymmetries
- Approx. Vigour Mid weight is from 0.9 to 2.5kg



Warranty Information

- 2 Years

Code	Display Name	Safe Working Load	Backrest Height	Backrest Width
PTB535032	Vigour Kidz - Mid Active Back Support For 12in W/C (12in High) For Mid Trunk Support	90 kg	12 in	11 in
PTB535048	Vigour - Mid Active Back Support For 14in W/C (14in High) For Mid Trunk Support	90 kg	14 in	13 in
PTB535056	Vigour - Mid Active Back Support For 15in W/C (14in High) For Mid Trunk Support	90 kg	14 in	14 in
PTB535064	Vigour - Mid Active Back Support For 16in W/C (16in High) For Mid Trunk Support	160 kg	16 in	15 in
PTB535072	Vigour - Mid Active Back Support For 17in W/C (17in High) For Mid Trunk Support	160 kg	17 in	16 in
PTB535080	Vigour - Mid Active Back Support For 18in W/C (17in High) For Mid Trunk Support	160 kg	17 in	17 in

Code	Display Name	Safe Working Load	Backrest Height	Backrest Width
PTB535096	Vigour - Mid Active Back Support For 20in W/C (17in High) For Mid Trunk Support	160 kg	17 in	19 in