

Liko ComfortVest Standing Sling (Mod. 95)



- Sit-to-stand sling designed especially for patients who must avoid underarm pressure
- Outer vest supports the patient behind the back and outside the arms, while the inner waist belt prevents gliding and keeps the patient from sliding out of the vest
- Self-regulating tensioning function
- Ideal choice for patients who are paralysed on one side or have other conditions that inhibit underarm pressure
- Available in Polyester

Warranty Information

- 1 Year

 Contract approved

Support @ Home

Function: Lifting persons

Category: Prescribed

Click [here](#) for more information

NDIS

Item Code:

05_123603166_0105_1_2

Hoist –

Composite:05_420_0105_1_2

Click [here](#) for more information

Compliance

AS/NZS ISO 10535:2011

Code	Display Name	Safe Working Load	Material
LSS394180	Comfort Vest Mod 95 Sling - Polyester - SMALL	200 kg	Polyester (Loops)
LSS394190	Comfort Vest Mod 95 Sling - Polyester - MEDIUM	200 kg	Polyester (Loops)
LSS394200	Comfort Vest Mod 95 Sling - Polyester - LARGE	200 kg	Polyester (Loops)

Code	Display Name	Safe Working Load	Material
LSS394210	Comfort Vest Mod 95 Sling - Polyester - XLARGE	200 kg	Polyester (Loops)