

Pulse Oximeter



- Economical, easy to use, and compact finger-tip pulse oximeter
- Simply attach the lightweight unit to the finger to obtain instant SPO2, pulse and Pleth Wave readings with the high visibility four-way OLED double colour display
- Integrated with SpO2 probe and processing display module
- Small in volume, light in weight and convenient in carrying



Contract approved

Support @ Home

Function: Not on S@H

Category: Not on S@H

Click [here](#) for more information

Warranty Information

- 1 Year

Code	Display Name
GMD332920	Pulse Oximeter - Finger Tip
GMD332925	Pulse Oximeter - Hand Held with Alarm