

Todd Adjustable Day Chair



- High back orthopaedic ergonomic comfort chair
- Built-in lumbar support, lateral support and knee bolster
- Fully adjustable seat depth and height
- Rubber stoppers for stability
- Easy to clean vinyl upholstery
- Powder coated steel frame
- Also available in pressure care model featuring pressure sensitive foams and 2-way stretch fabrics for the contact surfaces, combined with washable decorator fabrics on the outer surfaces
- Range of colours: Fabric: Tasman, Fern, Chocolate, Fawn, Rust. Vinyl: Sakara, Opal, Stonewash Blue, Stonewash Rush, Mist
- Specify colour at time of order

✓ Contracted DVA product

✓ Contract approved

ndis

Item Code:

05_188800051_0103_1_2

Click [here](#) for more information

Warranty Information

- 1 Year

Code	Display Name	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Safe Working Load
CHP193100	Todd Chair	395 - 510 mm	480 mm	750 mm	410 - 585 mm	130 kg