

TributeWrap Below Knee



- TributeWrap is an adjustable compression garment with foam technology that nestles into the tissue to promote skin health. It supports easy self-management of lymphoedema and oedema during low activity by comfortably adjusting to your patient's unique shape and lifestyle
- TributeWrap helps maintain gains made during therapy while allowing the flexibility to support further reduction. Simply slide it on, adjust the straps, and cover with Sleep Sleeve for comfortable evening compression

✔ Contracted DVA product

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

| Code | Display Name |
|-----------|-------------------------------------|
| FBS310030 | TributeWrap Below Knee Large Left |
| FBS310040 | TributeWrap Below Knee Large Right |
| FBS310050 | TributeWrap Below Knee Medium Left |
| FBS310060 | TributeWrap Below Knee Medium Right |
| FBS310070 | TributeWrap Below Knee Small Left |
| FBS310080 | TributeWrap Below Knee Small Right |