

Tandem Autosoft



The AutoSoft™ 90 infusion set is an all-in-one straight-set infusion system. The flexible 90-degree soft cannula comes pre-loaded in an automatic spring inserter, so one-handed insertions are quick and easy. It features:

- 90° insertion angle
- Soft Teflon cannula (6 or 9 mm length)
- Automatic insertion
- Reversible connector for quick reconnection
- Transparent window at insertion point
- Built-in needle cover
- Possible NDSS subsidy available.
- Fully funded for eligible DVA clients.

- ✓ Contracted DVA product
- ✓ Contract approved