


Sigvaris Ulcerx Compression Stocking Underliners



- Can be used in combination with the Sigvaris ULCER X Compression Stocking to treat severe venous conditions of the leg, including active ulceration
- Ideal for wound care compression therapy
- A delicate, slightly transparent stocking with a fine, elegant look
- Skin-friendly with optimal regulation of heat and moisture
- Relieves wound dressing discomfort during sleep or while recumbent
- Please specify size and colour at time of ordering

✔ Contract approved

Support @ Home

Function: Supporting and stimulating blood circulation

Category: Low risk

Click [here](#) for more information

Code	Display Name
PTS642400	Sigvaris - Ulcerx Underliners - Calf O/Toe (specify size & colour)